

ADJUSTABLE TORQUE WRENCH



Adjustable Torque Wrench # ATW

\$399

- ▶ Torque settings lock into select Location for accuracy and repeatability
- ▶ Ratchet clutch disengages (slips) when torque is applied to avoid over torquing

9 Torque Settings: 15, 20, 25, 30, 35, 40, 45, 50 & 60 N

Torque Control

Tightening with the correct torque minimizes the risk of the screw working loose and improves the precision during bridge fabrication. Torque Control can also reduce the risk for stress being built into the restoration and

Easy to use

Insert the screw or abutment into the implant with the appropriate driver and thumb knob manually and tighten with finger pressure. Final seat with the Precise Adjustable Torque Wrench with torque setting recommended by manufacturer. Pre-load the screw or abutment by using the torque wrench in the same way as a ratchet



DRIVER TIPS

Driver Insert Chart

IMPLANT SYSTEM	0.350 HCT	0.040 HSM	0.048 HNT	0.050 HTT	0.052 HS	0.062 HSL	0.080 HJ	0.080 HUL	0.080 HFD	STAR ITI
NobelBiocare	•		•		•	•			•	
Keystone Dental	•		•		•	•				
Steri-Oss				•						
Straumann/ITI										•
Friadent	•	•	•							
Impla-Med			•						•	
Implant Innovation /3i	•		•						•	
Imtec				•						
Implant Direct				•			•	•		
Zimmer				•						
BioHorizons	•			•						
Astra				•						
Camlog				•						
Innova				•						
PRICE	\$39	\$39	\$39	\$39	\$39	\$39	\$39	\$39	\$39	\$39

Abutment Torque Chart (Ncm)

Abutment	20 Ncm	30 - 32 Ncm	35 Ncm
Lifecore/Keystone	Lifecore/Keystone	Lifecore/Keystone	Lifecore/Keystone
3i/Biomet	3i/Biomet	3i/Biomet	3i/Biomet
Innova/Sybron	Innova/Sybron	Innova/Sybron	n/a
Astra	n/a	n/a	n/a
Friadent/Dentsply	Friadent/Dentsply	Friadent/Dentsply	n/a
Implant Direct	Implant Direct	Implant Direct	Implant Direct
n/a	Neoss	n/a	n/a
Zimmer	Zimmer	Zimmer	n/a
n/a	n/a	n/a	NobelBiocare
n/a	BioHorizons	n/a	n/a
n/a	n/a	n/a	Straumann/ITI
Imtec	n/a	n/a	n/a