

GERMAN-MADE MICRO SURGICAL LOUPES

LOWEST PRICES ANYWHERE (by 40 - 70%)

- Magnification enables SITTING or STANDING in a comfortable position, reducing eye, neck and back strain.

Our Quality EXCEEDS or BEATS that of our Competition!!!!

\$350.00
2.0 MAG only



| Mag | Working Distance | Depth of Field | Field of View | Weight | Pupillary Adjustmnt | Price |
|--------|------------------|----------------|---------------|------------|---------------------|-----------------|
| 2.0 S | 280mm | 125-203mm | 110mm | 68.5 grams | 51-76mm | |
| | 11" | 5-8" | 4.33" | 2.4 oz | 2-3" | |
| 2.0 R | 320mm | 125-203mm | 110mm | 68.5 grams | 51-76mm | |
| | 12.5" | 5-8" | 4.33" | 2.4 oz | 2-3" | \$350.00 |
| 2.0 L | 400mm | 125-203mm | 110mm | 68.5 grams | 51-76mm | |
| | 15.8" | 5-8" | 4.33" | 2.4 oz | 2-3" | |
| 2.0 XL | 460mm | 125-203mm | 110mm | 68.5 grams | 51-76mm | |
| | 18" | 5-8" | 4.33" | 2.4 oz | 2-3" | |
| Mag | Working Distance | Depth of Field | Field of View | Weight | Pupillary Adjustmnt | Price |
| 2.5 R | 350mm | 76-127mm | 76mm | 63.5 grams | 51-76mm | |
| | 13.8" | 3-5" | 3" | 2.2 oz | 2-3" | |
| 2.5 L | 480mm | 76-127mm | 76mm | 63.5 grams | 51-76mm | |
| | 19.0" | 3-5" | 3" | 2.2 oz | 2-3" | \$380.00 |
| Mag | Working Distance | Depth of Field | Field of View | Weight | Pupillary Adjustmnt | Price |
| 3.5 S | 350mm | 76-101mm | 65mm | 61.5 grams | 51-76mm | |
| | 13.8" | 3-4" | 2.5" | 2.1 oz | 2-3" | |
| 3.5 R | 420mm | 76-101mm | 65mm | 61.5 grams | 51-76mm | |
| | 16.53" | 3-4" | 2.5" | 2.1 oz | 2-3" | \$415.00 |
| 3.5 L | 500mm | 76-101mm | 65mm | 61.5 grams | 51-76mm | |
| | 19.68" | 3-4" | 2.5" | 2.1 oz | 2-3" | |
| Mag | Working Distance | Depth of Field | Field of View | Weight | Pupillary Adjustmnt | Price |
| 4.0 S | 350mm | 38-50mm | 57mm | 115 grams | 51-76mm | |
| | 13.8" | 1.5-2" | 2.25" | 3.9 oz | 2-3" | |
| 4.0 R | 460mm | 38-50mm | 57mm | 115 grams | 51-76mm | |
| | 18" | 1.5-2" | 2.25" | 3.9 oz | 2-3" | \$525.00 |

| Part # | ACCESSORIES | MSRP |
|-----------|--------------------------------|-------|
| LED | LED Power Pack * | \$670 |
| RM-STL- L | Safety Frames-Soft Touch | \$25 |
| RM-STL- S | Safety Frames-Soft Touch | \$25 |
| RM-TSF | Triwear Safety Frames | \$25 |
| RM-LC | Carrying Case | \$24 |
| RM-HAC | Hard Aluminum Case | \$35 |
| RM-HS | Head Straps | \$10 |
| RM-T-2 | Side Shields - Titanium(2/set) | \$6 |
| RM-L-2 | Side Shields - Aluminum(2/set) | \$6 |

** will clip onto MOST Optical Loupes

DVI / OROSCOPTIC USERS

SAVE 60 - 70%



LIFETIME Manufacturer Warranty on Optics

| Headlight | |
|-----------------------------|---------------------|
| Luminous intensity: | 140 lm |
| Illuminocity: | 70,000 lx (6500 fc) |
| Color temperature: | 6500K |
| Weight (w/o headband): | 14.5g (0.51 oz) |
| LED Lifetime: | 100,000 hours |
| Power Pack | |
| Internal battery voltage: | 4.8 V |
| Full battery charging time: | 120 min |
| Average operating time: | 5.5 hours |



\$670.00

Surgical

Solutions USA