

GERMAN-MADE MICRO SURGICAL LOUPES

LOWEST PRICES ANYWHERE (by 40 - 70%)

- Magnification enables **SITTING** or **STANDING** in a comfortable position, reducing eye, neck and back strain.

Our Quality **EXCEEDS** or **BEATS** that of our Competition!!!!

\$369.00
2.0 MAG only



Mag	Working Distance	Depth of Field	Field of View	Weight	Pupillary Adjustmnt	Price
2.0 S	280mm	125-203mm	110mm	68.5 grams	51-76mm	
	11"	5-8"	4.33"	2.4 oz	2-3"	
2.0 R	320mm	125-203mm	110mm	68.5 grams	51-76mm	
	12.5"	5-8"	4.33"	2.4 oz	2-3"	\$369.00
2.0 L	400mm	125-203mm	110mm	68.5 grams	51-76mm	
	15.8"	5-8"	4.33"	2.4 oz	2-3"	
2.0 XL	460mm	125-203mm	110mm	68.5 grams	51-76mm	
	18"	5-8"	4.33"	2.4 oz	2-3"	
Mag	Working Distance	Depth of Field	Field of View	Weight	Pupillary Adjustmnt	Price
2.5 R	350mm	76-127mm	76mm	63.5 grams	51-76mm	
	13.8"	3-5"	3"	2.2 oz	2-3"	
2.5 L	480mm	76-127mm	76mm	63.5 grams	51-76mm	
	19.0"	3-5"	3"	2.2 oz	2-3"	\$399.00
Mag	Working Distance	Depth of Field	Field of View	Weight	Pupillary Adjustmnt	Price
3.5 S	350mm	76-101mm	65mm	61.5 grams	51-76mm	
	13.8"	3-4"	2.5"	2.1 oz	2-3"	
3.5 R	420mm	76-101mm	65mm	61.5 grams	51-76mm	
	16.53"	3-4"	2.5"	2.1 oz	2-3"	\$429.00
3.5 L	500mm	76-101mm	65mm	61.5 grams	51-76mm	
	19.68"	3-4"	2.5"	2.1 oz	2-3"	
Mag	Working Distance	Depth of Field	Field of View	Weight	Pupillary Adjustmnt	Price
4.0 S	350mm	38-50mm	57mm	115 grams	51-76mm	
	13.8"	1.5-2"	3.25"	3.9 oz	2-3"	
4.0 R	460mm	38-50mm	57mm	115 grams	51-76mm	
	18"	1.5-2"	3.25"	3.9 oz	2-3"	\$529.00

Part #	ACCESSORIES	MSRP
LED	LED Power Pack *	\$649
RM-STL- L	Safety Frames-Soft Touch	\$25
RM-STL- S	Safety Frames-Soft Touch	\$25
RM-TSF	Triwear Safety Frames	\$25
RM-LC	Carrying Case	\$24
RM-HAC	Hard Aluminum Case	\$35
RM-HS	Head Straps	\$10
RM-T-2	Side Shields - Titanium(2/set)	\$6
RM-L-2	Side Shields - Aluminum(2/set)	\$6

** will clip onto MOST Optical Loupes

DVI / OROSCOPTIC USERS

SAVE 60 - 70%



LIFETIME Manufacturer Warranty on Optics

Headlight	
Luminous intensity:	140 lm
Illuminocity:	70,000 lx (6500 fc)
Color temperature:	6500K
Weight (w/o headband):	14.5g (0.51 oz)
LED Lifetime:	100,000 hours
Power Pack	
Internal battery voltage:	4.8 V
Full battery charging time:	120 min
Average operating time:	5.5 hours



\$649.00

Surgical

Solutions USA